



MC2000B[®]

Optical Chopper

SDK Manual





Table of Contents

Chapter 1	Introduction	1
Chapter 2	C++ Software Development Kit	2
2.1.	<i>cmd_library.h File Reference</i>	2
2.1.1.	Functions	2
2.1.2.	Function Documentation	4
Chapter 3	Python Software Development Kit	14
3.1.	<i>LK220_COMMAND_LIB Namespace Reference</i>	14
3.1.1.	Functions	14
3.1.2.	Function Documentation	15
Chapter 4	LabVIEW Software Development Kit	21
4.1.	<i>How to install</i>	21
4.2.	<i>How to find VI</i>	22
4.3.	<i>How to use</i>	22

Chapter 1 Introduction

User can start software development in C/C++ develop environment, LabVIEW and Python.

The software development interface can be found in the start menu.



or by clicking *Help* in software menu.

In this directory, you will find the support files for software development, as shown below.

is PC > Local Disk (C:) > Program Files (x86) > Thorlabs > MC2000B > Sample >		
Name	Date modified	Type
Thorlabs_MC2000B_C++ SDK	2/27/2020 11:20 PM	File folder
Thorlabs_MC2000B_LabVIEW SDK	2/27/2020 11:20 PM	File folder
Thorlabs_MC2000B_Python SDK	2/27/2020 11:20 PM	File folder
MC2000B SDK Manual	2/27/2020 10:47 PM	PDF File

Chapter 2 C++ Software Development Kit

User can start software development with **MC2000CommandLibWin32.dll** or **MC2000CommandLibWin64.dll** in C/C++ development environment which can be found in **Thorlabs_MC2000B_C++SDK** under \Document directory. The corresponding header file is also in **Thorlabs_MC2000B_C++SDK** under \Document directory.

Copy **MC2000CommandLibWin32.dll** or **MC2000CommandLibWin64.dll** to your program folder, and make sure the library file and exe file are in the same folder.

Below is the description of the header file **MC2000CommandLib.h**

2.1. cmd_library.h File Reference

2.1.1.Functions

- **MC2000COMMANDLIB_API int List** (char *sn)

list all the possible port on this computer.

- **MC2000COMMANDLIB_API int Open** (char *serialNo, int nBaud, int timeout)

open port function.

- **MC2000COMMANDLIB_API int IsOpen** (char *serialNo)

check opened status of port

- **MC2000COMMANDLIB_API int Close** (int hdl)

close current opened port

- **MC2000COMMANDLIB_API int Read** (int hdl, char *b, int limit)
- **MC2000COMMANDLIB_API int Write** (int hdl, char *b, int size)
- **MC2000COMMANDLIB_API int GetId** (int hdl, char *id)

Get the model number and firmware version.

- **MC2000COMMANDLIB_API int Set** (int hdl, char *c, int var)

- **MC2000COMMANDLIB_API int Get** (int hdl, char *c, char *d)

- **MC2000COMMANDLIB_API int SetTimeout** (int hdl, int time)

set time out value for read or write process.

- **MC2000COMMANDLIB_API int Purge** (int hdl, int flag)

Purge the RX and TX buffer on port.

- **MC2000COMMANDLIB_API int SetFrequency** (int hdl, int frequency)

Set the desired internal reference frequency.

- **MC2000COMMANDLIB_API int GetFrequency** (int hdl, int *frequency)

Get the internal reference frequency.

- **MC2000COMMANDLIB_API int SetBladeType** (int hdl, int type)

Set the blade type.

- **MC2000COMMANDLIB_API int GetBladeType** (int hdl, int *type)

Get the blade type.

- **MC2000COMMANDLIB_API int SetHarmonicMultiplier** (int hdl, int nharmonic)

Set Harmonic Multiplier applied to external reference frequency.

- **MC2000COMMANDLIB_API int GetHarmonicMultiplier** (int hdl, int *nharmonic)

Get Harmonic Multiplier applied to external reference frequency.

- **MC2000COMMANDLIB_API int SetHarmonicDivider** (int hdl, int dharmonic)

Set the Harmonic Divider applied to external reference frequency.

- **MC2000COMMANDLIB_API int GetHarmonicDivider** (int hdl, int *dharmonic)

Get the Harmonic Divider applied to external reference frequency.

- **MC2000COMMANDLIB_API int SetPhase** (int hdl, int phase)

Set the Phase adjust.

- **MC2000COMMANDLIB_API int GetPhase** (int hdl, int *phase)

Get the Phase adjust.

- **MC2000COMMANDLIB_API int SetEnable** (int hdl, int enable)

Set enable or disable state of device.

- **MC2000COMMANDLIB_API int GetEnable** (int hdl, int *enable)

Get enable or disable state of device.

- **MC2000COMMANDLIB_API int SetReference** (int hdl, int ref)

Set the reference mode.

- **MC2000COMMANDLIB_API int GetReference** (int hdl, int *ref)

Get the reference mode.

- **MC2000COMMANDLIB_API int SetReferenceOutput** (int hdl, int output)

Set the output reference mode.

- **MC2000COMMANDLIB_API int GetReferenceOutput** (int hdl, int *output)

Get the output reference mode.

- **MC2000COMMANDLIB_API int SetDisplayIntensity** (int hdl, int intensity)

Set the display intensity (1-10).

- **MC2000COMMANDLIB_API int GetDisplayIntensity** (int hdl, int *intensity)

Get the display intensity (1-10).

- **MC2000COMMANDLIB_API int GetReferenceFrequency** (int hdl, int *frequency)

Get the current supplied external reference frequency.

- **MC2000COMMANDLIB_API int Restore** (int hdl)

Restore the factory default parameters.

- **MC2000COMMANDLIB_API int GetVerbose** (int hdl, int *verbose)

Get the verbose mode.

- **MC2000COMMANDLIB_API int SetVerbose** (int hdl, int verbose)

Set the verbose mode.

- **MC2000COMMANDLIB_API int SetLanguage** (int hdl, int lang)

Set display language in the device.

- **MC2000COMMANDLIB_API int GetLocked** (int hdl, int *lock)

Get the lock state.

- **MC2000COMMANDLIB_API int GetCycleAdjust** (int hdl, int *cycle)

Get the cycle adjustment.

- **MC2000COMMANDLIB_API int SetCycleAdjust** (int hdl, int cycle)

Set the cycle adjustment.

- **MC2000COMMANDLIB_API int GetReferenceOutputFrequency** (int hdl, int * output)

Gets actual frequency blade spinning.

2.1.2.Function Documentation

MC2000COMMANDLIB_API int Close (int *hdl*)

Close current opened port

Parameters:

<i>hdl</i>	Handle of port.
------------	-----------------

Returns:

0: success; negative number: failed.

MC2000COMMANDLIB_API int Get (int *hdl*, char * *c*, char * *d*)

Set command to device according to protocol in manual and get the return string.

Make sure the port was opened successful before call this function.

Make sure this is the correct device by checking the ID string before call this function.

Parameters:

<i>hdl</i>	Handle of port.
<i>c</i>	input command string (<255)
<i>d</i>	output string (<255)

Returns:

0: success;
 0xEA: CMD_NOT_DEFINED;
 0xEB: time out;
 0xED: invalid string buffer;

MC2000COMMANDLIB_API int GetBladeType (int *hdl*, int * *type*)

Get the blade type.

Parameters:

<i>hdl</i>	Handle of port.
<i>type</i>	Pointer to the blade type.

Returns:

Negative number: failed.
 0: success; negative number: failed.
 -1: invalid string buffer;
 -2: time out;
 -3: time out;
 -4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetCycleAdjust (int *hdl*, int * *cycle*)

Get the cycle adjustment.

Parameters:

<i>hdl</i>	Handle of port.
<i>verbose</i>	Pointer to the cycle adjustment

Returns:

Negative number: failed.
 0: success; negative number: failed.
 -1: invalid string buffer;
 -2: time out;
 -3: time out;
 -4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetDisplayIntensity (int *hdl*, int * *intensity*)

Get the display intensity (1-10).

Parameters:

<i>hdl</i>	Handle of port.
<i>intensity</i>	Pointer to display intensity (1-10).

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetEnable (int *hdl*, int * *enable*)

Get enable or disable state of device.

Parameters:

<i>hdl</i>	Handle of port.
<i>enable</i>	Pointer to enable state 1: Enable 0: disable.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetFrequency (int *hdl*, int * *frequency*)

Get the internal reference frequency.

Parameters:

<i>hdl</i>	Handle of port.
------------	-----------------

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetHarmonicDivider (int *hdl*, int * *dharmonic*)

Get the Harmonic Divider applied to external reference frequency.

Parameters:

<i>hdl</i>	Handle of port.
<i>dharmonic</i>	Pointer to the harmonic divider applied to external reference frequency.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetHarmonicMultiplier (int *hdl*, int * *nharmonic*)

Get Harmonic Multiplier applied to external reference frequency.

Parameters:

<i>hdl</i>	Handle of port.
<i>nharmonic</i>	Pointer to the harmonic multiplier.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetId (int *hdl*, char * *id*)

Get the model number and firmware version.

Parameters:

<i>hdl</i>	Handle of port.
<i>id</i>	Pointer to id string buffer.

Returns:

0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetLocked (int *hdl*, int * *lock*)

Get the lock state.

Parameters:

<i>hdl</i>	Handle of port.
<i>verbose</i>	pointer to the lock state

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetPhase (int *hdl*, int * *phase*)

Get the Phase adjust.

Parameters:

<i>hdl</i>	Handle of port.
<i>phase</i>	pointer to the phase adjust.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetReference (int *hdl*, int * *ref*)

Get the reference mode.

Parameters:

<i>hdl</i>	Handle of port.
<i>ref</i>	Pointer to reference mode.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetReferenceFrequency (int *hdl*, int * *frequency*)

Get the current supplied external reference frequency.

Parameters:

<i>hdl</i>	Handle of port.
<i>frequency</i>	Pointer to the current supplied external reference frequency.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetReferenceOutputFrequency (int *hdl*, int * *freq*)

Get actual frequency blade spinning

Parameters:

<i>hdl</i>	Handle of port.
<i>freq</i>	Pointer to the actual frequency.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetReferenceOutput (int *hdl*, int * *output*)

Get the output reference mode.

Parameters:

<i>hdl</i>	Handle of port.
<i>output</i>	Pointer to the output mode.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;

-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int GetVerbose (int *hdl*, int * *verbose*)

Get the verbose mode.

Parameters:

<i>hdl</i>	Handle of port.
<i>verbose</i>	verbose mode

Returns:

Negative number: failed.
 0: success; negative number: failed.
 -1: invalid string buffer;
 -2: time out;
 -3: time out;
 -4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int IsOpen (char * *serialNo*)

Check opened status of port

Parameters:

<i>serialNo</i>	Serial number of the device to be checked.
-----------------	--

Returns:

0: port is not opened; 1: port is opened.

MC2000COMMANDLIB_API int List (char * *sn*)

List all the possible port on this computer.

Parameters:

<i>nPort</i>	port list returned string include serial number and device descriptor, separated by comma
--------------	---

MC2000COMMANDLIB_API int Open (char * *serialNo*, int *nBaud*, int *timeout*)

Open port function.

Parameters:

<i>serialNo</i>	Serial number of the device to be opened, use GetPorts function to get exist list first.
<i>nBaud</i>	bit per second of port
<i>timeout</i>	set timeout value in (s)

Returns:

Non-negative number: hdl number returned successfully; negative number: failed.

MC2000COMMANDLIB_API int Purge (int *hdl*, int *flag*)

Purge the RX and TX buffer on port.

Parameters:

<i>hdl</i>	Handle of port.
<i>flag</i>	

FT_PURGE_RX: 0x01

FT_PURGE_TX: 0x02

Returns:

0: success; negative number: failed.

MC2000COMMANDLIB_API int Read (int *hdl*, char * *b*, int *limit*)

Read string from device through opened port.

Make sure the port was opened successful before call this function.

Parameters:

<i>hdl</i>	Handle of port.
<i>b</i>	returned string buffer
<i>limit</i>	

ABS(limit): max length value of b buffer.

SIGN(limit) == 1 : wait RX event until time out value expired;

SIGN(limit) == -1: INFINITE wait event untill RX has data;

Returns:

Non-negative number: size of actual read data in byte; negative number: failed.

MC2000COMMANDLIB_API int Restore (int *hdl*)

Restore the factory default parameters.

Parameters:

<i>hdl</i>	Handle of port.
------------	-----------------

Returns:

Negative number: failed.

0: success; negative number: failed.

-1: invalid string buffer;

-2: time out;

-3: time out;

-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int Set (int *hdl*, char * *c*, int *var*)

Set command to device according to protocol in manual.

Make sure the port was opened successful before call this function.

Make sure this is the correct device by checking the ID string before call this function.

Parameters:

<i>hdl</i>	Handle of port.
<i>c</i>	input command string
<i>var</i>	length of input command string (<255)

Returns:

0: success;

0xEA: CMD_NOT_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

MC2000COMMANDLIB_API int SetBladeType (int *hdl*, int *type*)

Set the blade type.

Parameters:

<i>hdl</i>	Handle of port.
<i>type</i>	Desired blade type.

Returns:

Negative number: failed.

0: success; negative number: failed.

-1: invalid string buffer;

-2: time out;

-3: time out;
 -4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetCycleAdjust (int *hdl*, int *cycle*)

Set the cycle adjustment.

Parameters:

<i>hdl</i>	Handle of port.
<i>verbose</i>	the desired cycle adjustment

Returns:

Negative number: failed.
 0: success; negative number: failed.
 -1: invalid string buffer;
 -2: time out;
 -3: time out;
 -4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetDisplayIntensity (int *hdl*, int *intensity*)

Set the display intensity (1-10).

Parameters:

<i>hdl</i>	Handle of port.
<i>intensity</i>	Display intensity (1-10).

Returns:

Negative number: failed.
 0: success; negative number: failed.
 -1: invalid string buffer;
 -2: time out;
 -3: time out;
 -4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetEnable (int *hdl*, int *enable*)

Set enable or disable state of device.

Parameters:

<i>hdl</i>	Handle of port.
<i>enable</i>	Enable state 1: enable, 0: disable.

Returns:

Negative number: failed.
 0: success; negative number: failed.
 -1: invalid string buffer;
 -2: time out;
 -3: time out;
 -4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetFrequency (int *hdl*, int *frequency*)

Set the desired internal reference frequency.

Parameters:

<i>hdl</i>	Handle of port.
<i>id</i>	Desired internal reference frequency.

Returns:

Negative number: failed.
 0: success; negative number: failed.
 -1: invalid string buffer;

-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetHarmonicDivider (int *hdl*, int *dharmonic*)

Set the Harmonic Divider applied to external reference frequency.

Parameters:

<i>hdl</i>	Handle of port.
<i>dharmonic</i>	Desired harmonic divider applied to external reference frequency.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetHarmonicMultiplier (int *hdl*, int *nharmonic*)

Set Harmonic Multiplier applied to external reference frequency.

Parameters:

<i>hdl</i>	Handle of port.
<i>nharmonic</i>	Desired harmonic multiplier.

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetLanguage (int *hdl*, int *lang*)

Set display language in the device.

Parameters:

<i>hdl</i>	Handle of port.
<i>lang</i>	Language mode. 1: English, 0: Chinese

Returns:

Negative number: failed.
0: success; negative number: failed.
-1: invalid string buffer;
-2: time out;
-3: time out;
-4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetPhase (int *hdl*, int *phase*)

Set the Phase adjust.

Parameters:

<i>hdl</i>	Handle of port.
<i>phase</i>	Desired phase adjust.

Returns:

Negative number: failed.
0: success; negative number: failed.

- 1: invalid string buffer;
- 2: time out;
- 3: time out;
- 4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetReference (int *hdl*, int *ref*)

Set the reference mode.

Parameters:

<i>hdl</i>	Handle of port.
<i>ref</i>	Reference mode.

Returns:

- Negative number: failed.
- 0: success; negative number: failed.
- 1: invalid string buffer;
- 2: time out;
- 3: time out;
- 4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetReferenceOutput (int *hdl*, int *output*)

Set the output reference mode.

Parameters:

<i>hdl</i>	Handle of port.
<i>output</i>	Output mode.

Returns:

- Negative number: failed.
- 0: success; negative number: failed.
- 1: invalid string buffer;
- 2: time out;
- 3: time out;
- 4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int SetTimeout (int *hdl*, int *time*)

Set time out value for read or write process.

Parameters:

<i>hdl</i>	Handle of port.
<i>time</i>	time out value

MC2000COMMANDLIB_API int SetVerbose (int *hdl*, int *verbose*)

Set the verbose mode.

Parameters:

<i>hdl</i>	Handle of port.
<i>verbose</i>	Verbose mode. When verbose mode is set to 1, status messages are output on the USB

Returns:

- Negative number: failed.
- 0: success; negative number: failed.
- 1: invalid string buffer;
- 2: time out;
- 3: time out;
- 4: CMD_NOT_DEFINED;

MC2000COMMANDLIB_API int Write (int *hdl*, char * *b*, int *size*)

Write string to device through opened port.

Make sure the port was opened successful before call this function.

Parameters:






<i>hdl</i>	Handle of port.
<i>b</i>	input string
<i>size</i>	Size of string to be written.

Returns:

Non-negative number: number of bytes written; negative number: failed.

Chapter 3 Python Software Development Kit

Python 3.6 or above is required. User can import **MC2000B_COMMAND_LIB.py** to your python project, that's the wrapper for **MC2000B_COMMAND_LIB(MC2000BCommandLibraryWin32.dll** in C/C++ development environment). Copy **MC2000BCommandLibraryWin32.dll** to your program folder, and make sure the library file and **MC2000B_COMMAND_LIB.py** file are in the same folder. The “ **MC2000B_COMMAND_LIB_EXAMPLE.py**” is example code for how to use the python APIs.

	__pycache__	2/26/2020 1:28 PM	File folder	
	MC2000B_COMMAND_LIB.py	2/26/2020 2:04 PM	Python File	15 KB
	MC2000B_COMMAND_LIB_EXAMPLE.py	2/26/2020 2:13 PM	Python File	6 KB
	MC2000CommandLibWin32.dll	6/21/2016 10:01 AM	Application exten...	93 KB
	Thorlabs.MC2000B.CommandLibrary.Python.pyproj	2/25/2020 2:45 PM	Python Project	2 KB

User can also replace the reference win32 lib to x64 lib for 64-bit application.

3.1. LK220_COMMAND_LIB Namespace Reference

3.1.1.Functions

```

0)def MC2000BListDevices
1)def MC2000BOpen
2)def MC2000BIsOpen
3)def MC2000BClose
4)def MC2000BSetFrequency
5)def MC2000BGetFrequency
6)def MC2000BSetBladeType
7)def MC2000BGetBladeType
8)def MC2000BSetHarmonicMultiplier
9)def MC2000BGetHarmonicMultiplier
10)def MC2000BSetHarmonicDivider
11)def MC2000BGetHarmonicDivider
12)def MC2000BSetPhase
13)def MC2000BGetPhase
14)def MC2000BSetEnable
15)def MC2000BGetEnable
16)def MC2000BSetReference
17)def MC2000BGetReference
18)def MC2000BSetReferenceOutput
19)def MC2000BGetReferenceOutput
20)def MC2000BSetDisplayIntensity
21)def MC2000BGetDisplayIntensity
22)def MC2000BGetReferenceFrequency
23)def MC2000BRestore
24)def MC2000BGetVerbose
25)def MC2000BSetVerbose
26)def MC2000BSetLanguage
27)def MC2000BGetLocked
28)def MC2000BGetCycleAdjust
29)def MC2000BSetCycleAdjust
30)def MC2000BGetReferenceOutFrequency

```

3.1.2.Function Documentation

def MC2000B_COMMAND_LIB.MC2000BClose (hdl)

```
Close opened MC2000B device
Args:
    hdl: the handle of opened MC2000B device
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetBladeType (hdl, bladetype)

```
Get the blade type
Args:
    hdl: the handle of opened MC2000B device
    bladetype: pointer to the blade type
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetCycleAdjust (hdl, cycle)

```
Get the cycle adjustment.
Args:
    hdl: the handle of opened MC2000B device
    cycle: pointer to the cycle adjustment
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetDisplayIntensity (hdl, intensity)

```
Get the display intensity (1-10).
Args:
    hdl: the handle of opened MC2000B device
    intensity: pointer to display intensity (1-10)
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetEnable (hdl, enable)

```
Get enable or disable state of device
Args:
    hdl: the handle of opened MC2000B device
    enable: pointer to enable state 1:enable, 0:disable
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetFrequency (hdl, frequency)

```
Get the internal reference frequency
Args:
    hdl: the handle of opened MC2000B device
    frequency: get internal reference frequency
```

```
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetHarmonicDivider (hdl, dharmonic)

```
Get the Harmonic Divider applied to external reference frequency
Args:
    hdl: the handle of opened MC2000B device
    dharmonic: pointer to the harmonic divider applied to external reference frequency
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetHarmonicMultiplier (hdl, nharmonic)

```
Get Harmonic Multiplier applied to external reference frequency
Args:
    hdl: the handle of opened MC2000B device
    nharmonic: pointer to the harmonic multiplier
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetLocked (hdl, lock)

```
Get the lock state.
Args:
    hdl: the handle of opened MC2000B device
    lock: pointer to the lock state
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetPhase (hdl, phase)

```
Get the Phase adjust
Args:
    hdl: the handle of opened MC2000B device
    phase: pointer to the phase adjust
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetReference (hdl, ref)

```
Get the reference mode
Args:
    hdl: the handle of opened MC2000B device
    ref: pointer to reference mode 1:external, 0:internal
        reference mode high precision 0:InternalOuter 1:InternalInner 2:ExternalOuter
        3:ExternalInner
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetReferenceFrequency (hdl, frequency)

```
Get the current supplied external reference frequency.
Args:
    hdl: the handle of opened MC2000B device
    frequency: pointer to the current supplied external reference frequency.
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetReferenceOutFrequency (hdl, freq)

```
Get actual frequency of blade spinning
Args:
    hdl: the handle of opened MC2000B device
    freq: pointer to the actual frequency
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetReferenceOutput (hdl, output)

```
Get the output reference mode
Args:
    hdl: the handle of opened MC2000B device
    output: pointer to the output mode 1:actual, 0:target
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BGetVerbose (hdl, verbose)

```
Get the verbose mode.
Args:
    hdl: the handle of opened MC2000B device
    verbose: verbose mode
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BIsOpen (serialNo)

```
Check opened status of MC2000B device
Args:
    serialNo: serial number of MC2000B device
Returns:
    0: MC2000B device is not opened; 1: MC2000B device is opened.
```

def MC2000B_COMMAND_LIB.MC2000BListDevices ()

```
List all connected MC2000B devices
Returns:
    The MC2000B device list, each device item is [serialNumber, MC2000BType]
```

def MC2000B_COMMAND_LIB.MC2000BOpen (serialNo, nBaud, timeout)

```
Open MC2000B device
```

```
Args:
    serialNo: serial number of MC2000B device
    nBaud: bit per second of port
    timeout: set timeout value in (s)
Returns:
    non-negative number: hdl number returned Successful; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BRestore (hdl)

```
Restore the factory default parameters.
Args:
    hdl: the handle of opened MC2000B device
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetBladeType (hdl, bladetype)

```
Set the blade type
Args:
    hdl: the handle of opened MC2000B device
    bladetype: desired blade type
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetCycleAdjust (hdl, cycle)

```
Set the cycle adjustment.
Args:
    hdl: the handle of opened MC2000B device
    cycle: the desired cycle adjustment
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetDisplayIntensity (hdl, intensity)

```
Set the display intensity (1-10).
Args:
    hdl: the handle of opened MC2000B device
    intensity: display intensity (1-10)
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetEnable (hdl, enable)

```
Set enable or disable state of device
Args:
    hdl: the handle of opened MC2000B device
    enable: enable state 1:enable, 0:disable.
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetFrequency (hdl, frequency)

```
Set the desired internal reference frequency.
Args:
    hdl: the handle of opened MC2000B device
    frequency: desired internal reference frequency
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetHarmonicDivider (hdl, dharmonic)

```
Set the Harmonic Divider applied to external reference frequency
Args:
    hdl: the handle of opened MC2000B device
    dharmonic: desired harmonic divider applied to external reference frequency
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetHarmonicMultiplier (hdl, nharmonic)

```
Set Harmonic Multiplier applied to external reference frequency
Args:
    hdl: the handle of opened MC2000B device
    nharmonic: desired harmonic multiplier
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetLanguage (hdl, lang)

```
Set display language in the device.
Args:
    hdl: the handle of opened MC2000B device
    lang: language mode. 1: English, 0: Chinese
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetPhase (hdl, phase)

```
Set the Phase adjust
Args:
    hdl: the handle of opened MC2000B device
    phase: desired phase adjust
Returns:
    0: Success; negative number: failed.
```

def MC2000B_COMMAND_LIB.MC2000BSetReference (hdl, ref)

```
Set the reference mode
Args:
    hdl: the handle of opened MC2000B device
    ref: reference mode 1:external, 0:internal
    reference high pre
Returns:
```

```
0: Success; negative number: failed.
```

`def MC2000B_COMMAND_LIB.MC2000BSetReferenceOutput (hdl, output)`

```
Set the output reference mode
Args:
    hdl: the handle of opened MC2000B device
    output: output mode 1:actual, 0:target
Returns:
    0: Success; negative number: failed.
```

`def MC2000B_COMMAND_LIB.MC2000BSetVerbose (hdl, verbose)`


```
Set the verbose mode.
Args:
    hdl: the handle of opened MC2000B device
    verbose: verbose mode. When verbose mode is set to 1, status messages are output on the USB
Returns:
    0: Success; negative number: failed.
```

Chapter 4 LabVIEW Software Development Kit

The user can start software development with LabVIEW 2011 or later versions based on LabVIEW instrument driver mechanism. The supported files are in *LabVIEW SDK* under the **Sample** directory.

4.1. How to install


Unzip the zip file and **copy** to instr.lib folder under LabVIEW installation folder.



Destination folder: under %LabVIEW install path%\instr.lib


Typically, C:\Program Files (x86)\National Instruments\LabVIEW 2011\instr.lib

Note: LabVIEW 2011 or later LabVIEW versions are compatible.



4.2. How to find VI

VI Could be found under: Functions\Instrument I/O\Instrument Drivers\



4.3. How to use

1. From VI

Note: Before you open the SDK LabVIEW project, make sure the device has been connected to the computer.

Thorlabs_MC2000B.lvlib: Get Blade Type.vi (4815)







Handle in [11] — [3] Handle out
 error in (no error) [8] — [2] Blade Type
 — [0] error out

Blade Blades	Slots Dual Frequency Blades	Index	Single Frequency Blades	Higher Precision Inner
Outer	Inner		Outer	
MC1F2	2/100	0		4Hz - 200Hz
200Hz - 10KHz				
MC1F10	10	1	20Hz - 1KHz	
MC1F15	15	2	30Hz - 1.5KHz	
MC1F30	30	3	60Hz - 3KHz	
MC1F60	60	4	120Hz - 6KHz	
MC1F100	100	5	200Hz - 10KHz	

2. From VI tree


Some classic data flow in VI tree.

Use the Example Finder to find examples demonstrating the usage of this instrument driver. To launch Example Finder, select "Find Examples..." from the LabVIEW Help menu.

Open	Action/Status		Data		Close
					
		MC2000 GET REF FREQ	MC2000 SET LANG	MC2000 GET ID	MC2000 RESTOR
		MC2000 GET REF OUT FREQ			
		MC2000 GET LOCK			
				MC2000 GET BLADE	MC2000 SET BLADE
				MC2000 GET DISP	MC2000 SET DISP
				MC2000 GET FREQ	MC2000 SET FREQ
				MC2000 GET HARM DIWI	MC2000 SET HARM DIWI
				MC2000 GET HARM HULT	MC2000 SET HORM HULT
				MC2000 GET CYCLE ADJUST	MC2000 SET CYCLE ADJUST
				MC2000 GET ENABLE	MC2000 SET ENABLE
				MC2000 GET PHASE	MC2000 SET PHASE
				MC2000 GET OUTPUT	MC2000 SET OUTPUT
				MC2000 GET REF	MC2000 SET REF

3. From example

An examples show the classic usage. Example path: instr.lib\Thorlabs_MC2000B\Examples



Easy programming and detailed comment will help.



THORLABS
www.thorlabs.com
